

# 4s3™ Semi-Automatic Sheet Heat Sealer

Automatic sealing process for unmatched sealing consistency



## Fast - Simple - Flexible

- **Reproducible sealing quality and performance**  
Consistent and uniform sealing
- **Intuitive operation**  
Easy set-up of sealing temperature and time
- **Compatible with a wide range of plates and seals**  
Highest flexibility

# 4s3™ Semi-Automatic Sheet Heat Sealer

The 4s3™ is a semi-automatic sheet heat sealer which is compatible with a wide range of seals and plates of differing designs and heights.

With variable temperature and time settings, sealing conditions are easily optimised to produce a 100% tight seal, eliminating sample loss.

Plate and seal are placed on the holder, the "Operate" button is pressed and the drawer automatically closes. The sealing process is controlled by an electric mechanism providing consistent plate sealing, eliminating variation and giving reliable sealing every time.

- Automatic sealing process for unmatched sealing consistency
- Variable temperature and time settings for easy optimisation of sealing conditions
- Compatible with a wide variety of consumables
- Simple, single button operation
- Rapid heating element for fast start-up
- Real-time temperature display
- Auto stand-by function to save energy
- Auto switch off function for added safety
- RS-232 serial port enables full integration with robotic automation systems
- Full range of heat seals available
- Small footprint
- Low temperature sealing for use with adhesive seals



Real-time temperature display

## Features



### Flexibility of Plate Format

- Compatibility with a wide range of plate types with a maximum height of 60 mm
- The following plate formats can be used with the appropriate plate adapter:
  - Standard SBS footprint PP & PS plates & deepwell blocks
  - PCR plates (skirted, semi-skirted and non-skirted formats)



### Quick and Simple Installation

- Plug & play
- No set-up necessary



### Reproducible Sealing

- Reproducible sealing by applying a fixed pressure and accurate sealing time and temperature
- Adapters and sealing aids for seal optimization

### Adapters

- Deep well plates and most other skirted plates can be sealed on the standard plate support adapter (4ti-0615, included)
- For all other 96 well PCR plates, an additional plate support adapter (4ti-0625, optional) must be used
- Both adapters can also be combined for easier handling of shallow plates such as 384 well PCR plates
- For best sealing results with Random Access seals, 4 different adapters are available to offer the perfect plate support and seal alignment (see last page)




4ti-0615 Plate Support Adapter, Standard  
For skirted 96/384 well plates



4ti-0625 Plate Support Adapter, PCR 96  
For 96 well PCR plates

## Sealing Aids

- A weighted platen (4ti-0602) for use with deep well blocks can be ordered separately
- A sealing frame for use with all other plate designs is supplied when ordering the optional 4ti-0625 plate support adapter



4ti-0602  
Weighted Sealing Platen  
For use with deep well blocks



4ti-0612  
Sealing Frame  
For use with 4ti-0625



## Easy Handling

- The rapid heating element allows for a quick start-up within minutes
- The temperature is displayed in real-time so that the heating up process can be monitored
- Simple one button operation

## Operation in 6 Easy Steps



## RS-232 Port for Automation Integration

- The 4s3™ has an RS-232 port to enable it to be integrated into most automated systems
- SiLA compatible – quick and easy integration with other SiLA compliant devices using the SiLA driver

**SiLA** Rapid Integration®



## Energy Saving Functions

- When the unit is left idle it automatically switches into standby mode, during which the temperature of the heating element is reduced to 60°C to save energy
- When the sealer is left idle for longer periods it switches off automatically for added safety
- Standby times and switch-off modes can be adjusted to suit your needs

## Sealing Solutions

Heat sealing is the gold standard in plate sealing. It prevents sample loss and maximises sample security by ensuring a complete seal and preventing evaporation, leakage and contamination.

You have the option to choose your seal based on a wide variety of properties offered, including peelability, pierceability, gas permeability, optical clarity, temperature stability and solvent resistance.

Also see [www.4ti.co.uk/seal/heat-seals](http://www.4ti.co.uk/seal/heat-seals)

## FrameSeal™ Heat Sealing Frames

4titude® offers heat sealing sheets held within a disposable, rigid, stackable frame that are compatible with the 4s3™ heat sealer.

A FrameSeal™ can be placed either manually, or with a robotic gripper, onto a wide range of skirted and semi-skirted microplates (including PCR plates). After sealing, the frame can be removed by tearing along pre-cut lines to leave the heat seal in place on the plate surface.



## 4s3™ Heat Sealer & The Random Access Range

Random Access plates provide the greatest flexibility in PCR plates, affording handling of individual wells in a 96 well plate format.

The 4s3™ Heat Sealer allows for sealing of Random Access Plates in one step using Random Access heat seal sheets - 96 individual seal spots on a removable backing. This results in individually sealed tubes that are pierceable, allowing for easy sample access.

The corresponding Random Access adapter provides optimal support and alignment of Random Access seals for easy, reproducible sealing of individual wells.

For easy alignment, 4 pins are located on the adapter that the pinholes on the sealing sheets easily slide on to. Correct alignment can easily be checked thanks to the material between the seal spots being clear.



Random Access 96 Well Skirted PCR Plate (4ti-0960/RA)



Random Access Pierce Heat Seal Strong (4ti-05381/RA)




Place Random Access Plate on adapter



Align Random Access Seal with adapter



After sealing - remove seal backing liner



Individually sealed tubes



4ti-0613

Plate Support Adapter, RA, Low Profile

**4ti-0613**  
For low profile plates

**4ti-0613-1**  
For low profile Break-A-Way plates



4ti-0614

Plate Support Adapter, RA, Standard Profile

**4ti-0614**  
For standard profile plates

**4ti-0614-1**  
For standard profile Break-A-Way plates

### Ordering Information

| Code         | Description  |
|--------------|--|
| 4ti-0655     | <b>4s3™ Semi-Automatic Sheet Heat Sealer</b><br>Includes: Universal power supply, power cord, manual, standard plate support adapter (4ti-0615), 12 months parts and labour warranty |
| 4ti-0615     | 4s3™ Plate Support Adapter, Standard, supplied with the unit, for skirted 96 and 384 well plates   |
| 4ti-0616     | 4s3™ Plate Support Adapter, Roche 1536, for Roche 1536 well PCR plates   |
| 4ti-0625     | 4s3™ Plate Support Adapter, PCR 96 for 96 well PCR plates  |
| 4ti-0613     | 4s3™ Plate Support Adapter Random Access, Low Profile for low profile plates   |
| 4ti-0613-1   | 4s3™ Plate Support Adapter Random Access, Low Profile Break-A-Way for low profile FrameStar® Break-A-Way Plates  |
| 4ti-0614     | 4s3™ Plate Support Adapter Random Access, Standard Profile for standard profile plates   |
| 4ti-0614-1   | 4s3™ Plate Support Adapter Random Access, Standard Profile Break-A-Way for standard profile FrameStar® Break-A-Way Plates  |
| 4ti-0602     | 4s3™ Weighted Sealing Platen   |
| 4ti-0612     | 4s3™ Sealing Frame   |
| 4ti-0655-11  | Warranty extension, second year  |
| 4ti-0655-12  | Warranty extension, second and third year  |
| 4ti-0655-113 | Warranty extension, third year   |

### Technical Specifications

| Parameter                 | Specification      |
|---------------------------|--------------------|
| Dimension (W x L x H)     | 181 x 275 x 310 mm |
| Sealing Temperature Range | 60 – 200°C         |
| Sealing Time Range        | 0 – 10 sec         |
| Weight                    | 9 kg               |
| Power Supply              | 110/230 V          |
| Power Consumption         | 300 W (max)        |
| Power Saving Mode         | < 20 W             |
| Working Temperature Range | 15 – 30°C          |
| Operating Humidity (RH)   | 0 – 85             |
| Connection                | RS-232 serial port |

Distributed by: