


PCR Consumables

4titude offers a range of PCR consumables for low to medium throughput applications. Our tubes, tube strips and plates are moulded from virgin polypropylene under clean-room conditions. The same stringent QC procedures as for our FrameStar® range are used in production.



benefits


- Clean-room produced
- Thin-wall design
- Free of human DNA, RNases and DNases
- Leak tested



PCR Consumables

96well PCR Plate, Non-Skirted

- Universal fit for all 0.2ml thermal cyclers and most sequencers
- Standard tube height with 0.3ml working volume
- Black grid reference for easy sample identification
- Cut-off corner for plate orientation
- Clear wells for improved sample visibility
- Raised well rims to avoid cross contamination and to facilitate sealing
- Compatible with our domed and flat cap strips
- Heat sealing and adhesive sealing options available



CODE	DESCRIPTION	QTY
4ti-0750B	96well PCR plate, non-skirted	50 plates

96well PCR Plate, Skirted, Low Profile

- Low profile design, 0.2ml maximum working volume (when used with sealing sheets)
- Cut-off corner for plate orientation
- Alphanumeric grid reference
- Stackable
- Clear wells for improved sample visibility
- Raised well rims to avoid cross contamination and to facilitate sealing
- Compatible with our domed and flat cap strips
- Heat sealing and adhesive sealing options available
- Suitable for bar-coding



CODE	DESCRIPTION	QTY
4ti-0740C	96well PCR plate, skirted, low profile	50 plates

Strips of 8 PCR Tubes

- Suitable for all 0.2ml block thermal cyclers
- 0.25ml maximum recommended capacity
- Individually numbered tubes
- Can be cut into sections
- Available with domed and flat, optical caps



CODE	DESCRIPTION	QTY
4ti-0781	strips of 8 PCR tubes	125 strips
4ti-0780	strips of 8 PCR tubes and domed caps	125 strips
4ti-0784	strips of 8 PCR tubes and flat optical caps	125 strips

Strips of 8 Domed Caps

- Compatible with our tube strips and 96well PCR plates
- Easy to apply
- Large end tabs for easy removal
- Labelled for orientation



CODE	DESCRIPTION	QTY
4ti-0782	strips of 8 domed caps	125 strips
4ti-0752	strips of 8 domed caps	300 strips

Strips of 8 Flat Caps

- Compatible with our tube strips and 96well PCR plates
- Optically clear for fluorescence detection (e.g. Q-PCR)
- Easy to apply
- Large end tabs for easy removal
- Labelled for orientation



CODE	DESCRIPTION	QTY
4ti-0783	strips of 8 flat optical caps	125 strips
4ti-0751	strips of 8 flat optical caps	300 strips

PCR Consumables



Individual 0.2ml PCR Tubes

- Flat cap design
- Suitable for all 0.2ml block thermal cyclers
- Snap-shut cap
- 0.25ml maximum recommended capacity

CODE	DESCRIPTION	QTY
4ti-0790	0.2ml flat cap PCR tubes	1000 tubes



Related Products:

4titude® offers a range of 2 component design FrameStar® PCR plates for high throughput and low volume PCR as well as heat sealing and adhesive sealing options for all of our plates. Please contact 4titude or log on to www.4ti.co.uk for detailed information on these products.



Unit 4bc, Jayes Park Courtyard,
Ockley, Surrey RH5 5RR, UK
Tel: +44 (0)1306 621111
Fax: +44 (0)1306 621162
Web: www.4ti.co.uk
Email: sales@4ti.co.uk